


· Für HF- und RF-Generatoren geeignet  
*Applicable for HF- and RF-generators*

· Hochwertige Isolation  
*Valuable insulation*

· Europäischer Flachstecker / Führungsstift  
*European flat plug / guiding pin*



VPE / unit: 1  
Umax: 500 Vp



| Spitzenbreite   tip width | Nadel / Needle | 0.5 mm       | 1.0 mm       | 2.0 mm       | 2.0 mm Karo  |
|---------------------------|----------------|--------------|--------------|--------------|--------------|
|                           |                |              |              |              |              |
| Artikel-Nr.   item code   | -              | -            | <b>95042</b> | -            | -            |
| L1                        | -              | -            | 145 mm       | -            | -            |
| L2                        | -              | -            | 8 mm         | -            | -            |
| Artikel-Nr.   item code   | <b>95036</b>   | <b>95013</b> | <b>95014</b> | -            | -            |
| L1                        | 165 mm         | 165 mm       | 165 mm       | -            | -            |
| L2                        | 6 mm           | 8 mm         | 8 mm         | -            | -            |
| Artikel-Nr.   item code   | <b>95187</b>   | <b>95089</b> | <b>95090</b> | <b>95088</b> | -            |
| L1                        | 180 mm         | 180 mm       | 180 mm       | 180 mm       | -            |
| L2                        | 6 mm           | 8 mm         | 8 mm         | 8 mm         | -            |
| Artikel-Nr.   item code   | <b>95103</b>   | <b>95104</b> | <b>95001</b> | <b>95007</b> | <b>95058</b> |
| L1                        | 195 mm         | 195 mm       | 195 mm       | 195 mm       | 195 mm       |
| L2                        | 8 mm           | 8 mm         | 8 mm         | 8 mm         | 8 mm         |
| Artikel-Nr.   item code   | <b>95128</b>   | <b>95122</b> | <b>95024</b> | <b>95125</b> | -            |
| L1                        | 225 mm         | 225 mm       | 225 mm       | 225 mm       | -            |
| L2                        | 8 mm           | 8 mm         | 8 mm         | 8 mm         | -            |